
# **KaVo IXS™ Holders**

User Manual

## **Supports KaVo IXS**

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
# KaVo IXS™ Holders

## User Manual


# Chapter 1

## KaVo IXS Intraoral Positioning Device

The KaVo IXS Positioning Device is specifically designed to support the sensor and align it with the X-ray beam and with the teeth of the upper and lower jaws.


## IXS Sensor System Assembly




### Anterior

Apply a new disposable sheath prior to use


### Posterior


### Bitewing

#### Horizontal

#### Vertical

---

## Positioning the IXS Sensor

1. Assemble the rings, bars and holders so that the entire sensor area is visible through the ring. This applies to horizontal or vertical Bitewing holders, Periapical (anterior or posterior) holders, and hand-held holders.

2. Place the sheath over the sensor.

3. Attach the sensor to the positioner. The sensor will only fit into the holder one way.
4. Place the cable into the groove on the holder whenever possible. This protects it when the patient closes.

The sensor may alternatively be positioned by hand and held in place by the patient with one finger just like conventional X-ray film. If the positioning device is not used, a wad of cotton wool may be useful for aligning the sensor parallel to the tooth. It may be best to position the sensor by hand in the small mouths of children.

---

### **WARNING**

Under no circumstances should the dental professional hold the sensor by hand during the X-ray exposure.

---

Note however, that image quality might be inferior with respect to one attained using an aiming device.


---

### **CAUTION**

Biting down on the sensor or cable will damage the housing and/or cable and lead to a sensor malfunction.

---

5. Position the X-ray generator as usual for X-ray film.  
Use of the paralleling technique is highly recommended.

**Important:** Verify connectivity of the KaVo IXS to the computer along with having a green driver icon before attempting to take an X-ray.

**CAUTION** \_\_\_\_\_

For maximum image quality, make sure the patient remains still during the exposure.

---


**WARNING** \_\_\_\_\_

Take the necessary steps to protect yourself from radiation. For proper operator positioning, refer to the 'Instructions for Use' of your intraoral X-ray equipment.

---

6. Acquire the image by pressing the X-ray button on your radiography system.

**Important:** It is important that the X-ray beam hits the whole sensor surface.

7. Carefully remove the sensor from the oral cavity following the exposure.
8. Grasp the sensor and positioning device and remove the sensor from under the retaining clips.
9. If needed, grasp the USB connector and pull it straight out of the port on the computer.

---

## Pedodontic Placement Techniques

When you are placing a sensor in the mouth of a child, you will often face the same challenges that you face with placing film in their mouths. Keep in mind that although the sensor is larger than a film pack, it is also smoother. This will work to your advantage.

Here are some techniques that will help you succeed in using the sensor in children when standard techniques fail. You may be familiar with many of them, since they work not only for the KaVo IXS size 1 sensor, but for film as well.

To begin, you should set the child up for X-rays both physically and psychologically; that is, ready their mouths and their minds. Use diversion when necessary. In addition, consider the size of the pedo film relative to the sensor. Also consider what information you need to gather from the X-ray and use that part of the sensor that will enable you to capture the information you need.

- Avoid saying the word *tongue*, since it might cause the tongue to act involuntarily.
- Have the child focus on the end result — an instant image of their teeth. If they can see the monitor, this will give them incentive to cooperate. Most children are very comfortable with computers and will be impressed to see an X-ray of their teeth on the screen.
- Place your finger in the areas where you plan on placing the sensor.
- Let the patient hold the sensor in the mouth (without biting it) to get used to the feel of it.
- Physically divert focus to something other than the sensor. Have the child place his or her finger on their arm, hold the holder, or make a certain sound.

## Cable Protection

Depending on the holder type, there are cable grooves to help guide and protect the sensor cable out of the oral cavity. Place the sensor cable into this groove and use the Aiming Bar Cable Clips for protection. This works well when the sensor is placed into the holder and then the entire assembly is sleeved.

On those occasions where it is necessary to sleeve the sensor before placing in the holder, due to the thickness of the barrier, the cable may not stay in the groove. In this case, verify that the cable is adjacent to the holder, either to the side that is to the anterior or posterior of the holder. Since the holder will not allow for full closure, the cable is protected. Again, please verify that the cable is not between the holder and the teeth as it may become crimped upon the patient closing.

You must always ask patients to close **lightly** so that they have only the slightest of contact. Do not allow patients to bite down on the cable. If you feel that you will not have the patient's full cooperation, place an object such as a folded gauze or cotton roll on the tab, to hold the teeth open just enough to keep the cable from harm. You may also place the cable in an area where there are missing teeth or in an area where, anatomically, there is an open space or lack full occlusion.

---

## Recommended Disinfecting Solutions

### CAUTION

---

Holders should be disinfected in an autoclave or with a cleaning agent after each patient.

---

Decontaminate the sensor, cable and holders in accordance with CDC – or your country's (e.g., OSAP) standards – for infection control. KaVo recommends wiping the holders with one of the following cleaning agents.

- CaviCide® (Manufactured by Metrex)
- CaViWipes™ (Distributed by Kerr)
- ProSpray™ (manufactured by Certol)
- Sani-Cloth® Plus, HB (manufactured by PDI)

### DO NOT:

- Clean the sensor using non appropriate tools.

The following solutions can be used to decontaminate the holders:

The sensor holders can be disinfected in a steam autoclave using distilled water at temperature of 275°F (125°C) at 216 kPa for 12 minutes when bagged.

However, certain restrictions apply:

- Run the appropriate autoclave cycle according to the manufacturer's specifications for your specific unit.
- Use the plastic or rubber modes, if available.
- Always avoid direct contact of the holders with metal trays, instruments, and heating elements.
- Always place holders in autoclave bags.
- Always situate bagged holders into the autoclave so that they are furthest from the heat source.
- Do not use a chemical autoclave.


Exposing the holders to hot metal and placing them close to heating elements will reduce the lifetime of the holders.

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**Dental Imaging Technologies Corporation**
**Kaltenbach & Voigt GmbH**