A-dec® 300



The A-dec Difference Built for Life

At A-dec, no detail is too small, right down to the bolts that go into our products. That's why we manufacture a large majority of our finished dental equipment from raw materials: an essential factor in ensuring quality. Because while these parts may not be visible to the eye, they make the difference in the life, comfort and functionality of each and every product that bears our name. And to your business, that's everything.







Control block – Patented design promotes clean treatment water with every activation.

Lift cylinder – A-dec built and engineered with robust hydraulics and fewer parts.

Chair position sensor – A-dec technology uses precise electronic sensors and digital signals.





A-dec 300: Evolved

Design

A solid baseplate is the foundation for a sturdy chair. Along with sure footing, the A-dec 300 quietly and smoothly positions patients with a hydraulic lift cylinder that's engineered for a 20 year life.

Value

Designed as a modular solution, A-dec 300 is a whole system that lets you pick and choose the features you want, within the price point you need.

Flexibility

Traditional or Continental, the choice is yours. Whether you select a dedicated left/right delivery or ambidextrous Radius[®] configuration, handpieces and clinical devices are always within natural reach.

Comfort

Pressure mapped, contoured surfaces and virtual pivot back ensure a relaxed and supported patient experience. Comfortably accommodates patients up to 400 pounds (181 kg.).

Access

With your knees tucked under the ultra-thin and flexible backrest, you have optimal access to the oral cavity, whether you work at 9 o'clock or even 1 o'clock. A dual-articulating headrest and ample shoulder support keeps patients resting comfortably.



Ergonomic and Economic

Who says ultimate comfort isn't attainable? Every patient chair, and dentist or assistant stool we design, is backed by decades of research, dialogue, tests and improvements. We never stop looking at every simple curve, smooth line and inherent contour, to make sure it delivers ultimate efficiency and comfort. Rest assured: The affordable A-dec 300 is no exception.

Works the Way You Do

Radius or post mount, you choose. Both options place instrumentation within reach for uninterrupted patient care.



A-dec 300 Radius System

Switch between right or left-handed operators with ease. Single- or dual-Radius mounting rotates around the front of the chair for better patient access during procedures.



Lever Foot Control

A-dec's lever foot control option keeps you in control. Intuitive wet/dry sideto-side operation and chip blower/ accessory button lets you focus on patient treatment without looking away. With minimal foot movement, you get precise modulation of both electric and pneumatic handpieces.



A-dec 300 Radius System with Manual Height Adjustment

The A-dec 336 Traditional delivery system allows you to add an optional saliva ejector. Pair it with an A-dec 351 or A-dec 551 assistant's arm and create a second saliva ejector for the operatory.



Saliva Ejector Add to the A-dec 336 Traditional delivery system for an all-in-one hygiene package.

So Many Choices



Sewn upholstery with armrests



Formed upholstery without armrests

Upholstery and Armrests

It all starts with a great chair. You have plenty of options, from plush sewn or easy-to-clean formed upholstery, to two-position armrests or no armrests.



Traditional



Continental

Delivery Systems

How do you and your staff work? Traditional and Continental delivery systems can integrate a quad-volt intraoral light source and two ancillaries.





Standard touchpad

Disc foot control



Deluxe touchpad

Touchpad

Fingertip access simplifies chair adjustments and provides quick access to dental light and cuspidor control. Optional deluxe touchpad adds controls for integrated clinical devices such electric motors, cameras, and scalers.



Lever foot control

Foot Controls

Disc style with chip blower optionor A-dec lever style for precise modulation of electric and pneumatic handpieces.



Fixed glass cuspidor



Rotating vitreous china cuspidor

Cuspidors

What's your preference, glass or vitreous china? Fixed or rotating? Both bowls include timed cup fill and bowl rinse functions.







Assistant's Holders

Place the assistant's clinical devices right where they are needed most. Adjustable 3- or 4-position holders are available as a chair mount, cuspidor mount or telescoping arm.



A-dec 300 LED



Radius monitor mount



Light post monitor mount



Light post monitor mount with link arm

Monitor Mounts

Consult. Educate. Review. Monitor mounts pivot to accommodate viewing from the seated and supine positions. Tailor your options, including locking style or floating pivot, positioning handle, and multiple mount locations.

A-dec 500 LED

Dental Lights

Accurate lighting is essential. A-dec LED lights provide balanced intensity and a uniform light pattern without shadows. Automatic on/off when chair moves to entry or exit position.



Color makes all the difference in creating an environment that reflects you. That's why A-dec offers a wide variety of colors and selected options to coordinate formed and sewn upholsteries.

Sewn Upholstery



Actual colors may vary. See your authorized A-dec Dealer for color samples and the most current product information.



Order Your Samples at a-dec.com/InspireMe

Specifications

A-dec 311 Dental Chair

| Range of lift height | 13.75" to 31.5" (349 mm to 800 mm); Low point 14.75" (375 mm) when dual mounting | |
|----------------------------------|--|--|
| Thin flexible backrest | 1.25" (32 mm), formed | |
| Chair swivel | 30° each side of center | |
| Drive system | Hybrid drive (electromechanical tilt, hydraulic base) | |
| Chair back and base stop switch | Yes, stops downward movement when obstruction is detected | |
| Headrest | Dual-articulating and gliding (knob style) | |
| Integrated 300-watt power supply | Yes | |
| Foot switch | Optional | |
| Mounting capability | Floating - 2" post mount or Radius (single or dual mounting) | |
| Upholstery options | Sewn or formed | |
| | | |

Delivery Systems

| | A-dec 300 Continental | A-dec 300 Traditional |
|-------------------------------------|--|--|
| Control block | Standard 3-position, optional 4-position | Standard 3-position, optional 4-position |
| Left/right conversion | Yes | Yes |
| Tubing options | Vinyl or silicone | Vinyl or silicone |
| Handpiece holder positions | 4 | 4 |
| Brake handle(s) | 1 standard, 2nd optional | 1 standard, 2nd optional |
| Touchpad options | None, standard or deluxe | None, standard or deluxe |
| Quad-voltage intraoral light source | Optional | Optional |
| Foot control | Standard, chip blower or lever | Standard, chip blower or lever |
| Arm system | Balanced flexarm with air brake A-dec 333 - Radius mount A-dec 335 - Support mount | Balanced flexarm with air brake A-dec 332 - Radius mount A-dec 334 - Support mount Manual height adjustment A-dec 336 - Radius mount |
| Integrated saliva ejector | N/A | Optional on A-dec 336 |











A-dec Headquarters

2601 Crestview Drive Newberg, Oregon 97132 USA Tel: 1.800.547.1883 within USA/CAN Tel: +1.503.538.7478 outside USA/CAN Fax: +1.503.538.0276 a-dec.com

A-dec Australia

Tel: 1.800.225.010 within AUS Tel: +61.(0)2.8332.4000 outside AUS

A-dec China

Tel: 400.600.5434 within CN Tel: +86.571.89026088 outside CN

A-dec United Kingdom Tel: 0800.233285 within UK Tel: +44.(0).24.7635.0901 outside UK The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

KaVo is a trademark or registered trademark of KaVo Dental GmbH in the US and other countries.

©2021 A-dec Inc. All rights reserved. 85.6152.00/BB/1.5M/6-21/Rev M