

Planmeca PlanMill® 50 S

The versatile Planmeca PlanMill 50 S unit is a wet or dry mill capable of handling discs and blocks. Equipped with a high-speed spindle and an automatic changer for 12 tools, the powerful milling unit has been designed specifically for high-quality restorations and prosthetics.



Product features

- **5-axis** milling unit
- Ability to use **discs or blocks**, milling and grinding
- Choice of **wet or dry** processing
- Capable of milling a single crown to full arch restoration
- **Broad range of materials** available including **zirconia**, IPS e.max®, glass ceramics, temporary materials and much more
- **12 tool** automatic changer
- Measurement and tool check after each processing step
- Integrated PC with touchscreen
- Easy user maintenance



PLANMECA
Protection Program

*Keep your investments in peak performance
while you deliver exceptional dentistry.*



IPS e.max® is a registered trademark of Ivoclar Vivadent AG

PLANMECA