


Planmeca ProSensor[®] HD

Usability · Integration · Connectivity

Usability is standard with our Planmeca ProSensor HD. It can be integrated with our Planmeca ProX[™] intraoral X-ray unit or connected via Ethernet, a Plug & Play USB port, or with 3rd-party vendors open architecture systems. The elegant control box features a distinguished design equipped with a 4-color LED light providing instant visual feedback of the imaging status. Magnetically connected sensors are designed for easy one-handed use and portability, while the sensor's hermetically sealed design offers simplified infection control. All combine to help you work quickly, easily and efficiently.

Planmeca ProSensor HD advantages

- Available in three sizes (0, 1, 2), slim with rounded corners
- Magnetic sensor connection for one-handed use and unsurpassed damage protection
- Feature a theoretical resolution of 33 lp/mm
- Color coded LED indicator lights to alert the operator for sensor connectivity
- Durable Kevlar construction and bite-resistant exterior for a lasting investment
- 55lb. pull force helps maintain connection to unit
- Can be fully integrated into Planmeca ProX[™] with all wiring and connections concealed - no exposed connections
- Third-party vendor connectivity



What We Offer



Pairs perfectly with
Planmeca ProX[®]



High definition
image quality



Hermetically sealed for
easy infection control

Planmeca ProX™

Quality Engineered...Diagnostically Focused



Pairs perfectly with
Planmeca ProSensor® HD



0.4mm Focal spot for
exceptional image quality



Integrated control box for
easy accessibility

The Planmeca ProX intraoral x-ray system adapts to both short-cone and long-cone imaging techniques. Versatility is a Planmeca ProX standout feature with five different arm lengths: 60", 65", 70", 75", 80", and variable kV/mA to meet your imaging needs. Additionally, Planmeca ProX offers a smart control for maintaining consistency of radiographs whenever imaging conditions change, and a self-diagnostic control system that monitors all functions and displays error messages to assist in the correct use and service of the unit.

Planmeca ProX features

- Quick and easy to use: pre-programmed quick settings
- Integrates with Planmeca ProSensor® HD Intraoral Sensor
- Easy integration/workflow with Planmeca Romexis® open architecture software
- **Wide-range of installation possibilities**



For a complimentary in-office consultation, please call
Planmeca at (800) 718 4598 or visit us on the web at www.planmeca.com

PLANMECA

© Planmeca U.S.A. Inc. 2020 All Rights Reserved